


LAB

Identification of *Campylobacter jejuni* and *C. coli*, *Enterococcus faecium* and *E. faecalis* by PCR

Professor Rungtip Chuanchuen DVM MS PhD

Faculty of Veterinary Science, Chulalongkorn University

Why PCR for *Campylobacter* and *Enterococcus*?



Biochemical and phenotypic tests do not usually produce correct species differentiation for all *Campylobacter* and *Enterococcus* isolates. Molecular testing (e.g. PCR techniques) should be applied to confirm the species identity.

Amplification of *ddl* (D-alanine-D-alanine ligase) gene (**Dutka-Malen et al. (1995)**)

Target	Primer	PCR product (bp)
<i>E. faecalis</i>	FW 5'-ATCAAGTACAGTTAGTCTT-3'	941
	RV 5'-ACGATTCAAAGCTAACTG-3'	
<i>E. faecium</i>	FW 5'-GCAAGGCTTAGAGA-3'	550
	RV 5'-CATCGTGTAAGCTAACTTC-3'	


1. Pre-denaturation	94°C	2 min
2. Denaturation	94°C	1 min
3. Annealing	54°C	1 min
4. Extension	72°C	1 min
Repeat 2 – 4 for 30 cycles		
5. Final extension	72°C	10 min
6. Storage	4°C	∞

Amplification of **sodA** (superoxide dismutase) gene (**Jackson et al. (2004)**)

Target	Primer	PCR product (bp)
<i>E. faecalis</i>	FW 5'-ACTTATGTGACTAACTTAACC-3'	360
	RV 5'-TAATGGTGAATCTTGGTTTGG-3'	
<i>E. faecium</i>	FW 5'-GAAAAACAATAGAAGAATTAT-3'	215
	RV 5'-TGCTTTTTTTGAATTCTTCTTTA-3'	

1. Pre-denaturation	95°C	4 min
2. Denaturation	95°C	30 sec
3. Annealing	55°C	1 min
4. Extension	72°C	1 min
Repeat 2 – 4 for 30 cycles		
5. Final extension	72°C	7 min
6. Storage	4°C	∞

Identification of *Enterococcus* species using multiplex PCR (Jackson et al., 2004)



Lanes 4 to 7: L4 = *E. faecalis* (360 bp); L5 = *E. durans* (295 bp);
L6 = *E. faecium* (215 bp); L7 = *E. malodoratus* (134 bp)


>, enterococcus genus-specific amplicon based on 16S rRNA primers (733 bp)

Multiplex PCR for
Campylobacter genus,
C. jejuni, and *C. coli*
(Wang et al., 2002)

Target	Primer	PCR product (bp)
<i>Campylobacter</i> spp. (23S rRNA)	FW 5'-TATACCGGTAAGGAGTGCTGGAG-3' RV 5'-ATCAATTAACCTTCGAGCACCG-3'	650
<i>C. jejuni</i> (<i>hipO</i>)	FW 5'-ACTTCTTTATTGCTTGCTGC-3' RV 5'-GCCACAACAAGTAAAGAAGC-3'	323
<i>C. coli</i> (<i>glyA</i>)	FW 5'-GTAAAACCAAAGCTTATCGTG-3' RV 5'-TCCAGCAATGTGTGCAATG-3'	126

1. Pre-denaturation	95°C	6 min
2. Denaturation	95°C	30 sec
3. Annealing	57°C	30 sec
4. Extension	72°C	30 sec
Repeat 2 – 4 for 30 cycles		
5. Final extension	c	7 min
6. Storage	4°C	∞

PCR for identification of *Campylobacter* species



Conventional PCR



Carry out PCR under
sterile conditions.

Flow of a conventional PCR experiment



1

Make a plan




Prepare PCR components



Set PCR reaction



Run thermocycler




Run electrophoresis



PCR product visualization

1



Identification of *Campylobacter* species



1. Sample C1
2. Sample C2

2

Identification of *Enterococcus* species





1. Sample E1
2. Sample E2

Make a plan for your PCR


Materials

- Molecular marker (Ladder 100bp)
- Electrophoresis buffer (TAE or TBE)
- Tips (filter) for pipettes 1 μ L to 1000 μ L
- Agarose
- PCR components
- TE buffer
- Tris HCl buffer
- Crushed ice
- Mineral oil (if necessary)
- etc



Equipment

- PCR thermocycler
- Pipettes for 1 μ L to 1000 μ L
- Eppendorf tubes
- PCR tubes
- Electrophoresis unit
- Microwave
- Autoclave
- Eppendorf centrifuge (PCR tubes)
- Photo camera
- UV-transilluminator
- Water bath



etc

Flow of a conventional PCR experiment



1

Make a plan



2

Prepare PCR components




3

Set PCR reaction



4

Run thermocycler



5


Run electrophoresis



6

PCR product visualization




Components of a PCR reaction



- Template DNA
- Thermostable DNA polymerase
- Primers
- dNTPs
- Buffer and salts (KCl, MgCl₂)
- PCR grade water
- Optional: DMSO, Formamide

Template DNA preparation

Whole cell boiling methods



✓ Suspend a loopful of culture in 100 μ l of sterile DW.

✓ Boiled 10 min at 100°C

✓ Place on ice



✓ Centrifuged 5 min at 12,000-13,000xg

✓ Store at -20°C

✓ Transfer supernatant to a new tube



Flow of a conventional PCR experiment



1

Make a plan



2

Prepare PCR components




3

Set PCR reaction



4

Run thermocycler



5

Run electrophoresis



6

PCR product visualization

Setting PCR reaction



Keep and thaw all reagents on ice.

Assemble reaction mix into 50 μ L volume in a thin walled 0.2 mL PCR tubes.

Add reagents in following order:

Water \rightarrow buffer \rightarrow dNTPs \rightarrow MgCl_2 , \rightarrow template DNA \rightarrow primers \rightarrow Taq polymerase.




Gently mix by tapping tube. Briefly centrifuge to settle tube contents.

Setting PCR reaction

1

Enterococcus species



Sample E1



Sample E2



Faecium



Faecalis

Positive control



Negative control

2

Campylobacter species



Sample C1



Sample C2



C. jejuni




C. coli

Positive control



Negative control

Setting PCR reaction





Primers for multiplex PCR

Target	Primer	PCR product (bp)
<i>C. jejuni hipO</i>	<u>FW</u> 5'-ACTTCTTTATTGCTTGCTGC -3'	323
	<u>RV</u> 5'-GCCACAACAAGTAAAGAAGC-3'	
<i>C. coli glyA</i>	<u>FW</u> 5'-GTAAAACCAAAGCTTATCGTG-3'	126
	<u>RV</u> 5'-TCCAGCAATGTGTGCAATG-3'	
<i>C. jejuni 23s rRNA</i>	<u>FW</u> 5'-TATACCGGTAAGGAGTGCTGGAG -3'	650
	<u>FW</u> 5'-ATCAATTAACCTTCGAGCACCG -3'	

Wang et al. (2002)



Primers for multiplex PCR

Target	Primer	PCR product (bp)
<i>E. faecalis</i>	<u>FW</u> 5'-ACTTATGTGACTAACTTAACC -3'	360
	<u>RV</u> 5'-TAATGGTGAATCTTGGTTTGG-3'	
<i>E. faecium</i>	<u>FW</u> 5'-GAAAAACAATAGAAGAATTAT-3'	215
	<u>RV</u> 5'-TGCTTTTTTGAATTCTTCTTTA-3'	

Jackson et al. (2004)

1

Enterococcus spp.

Component	Stock conc.	Volume in 1 rxn (μL)	No. of rxn	Volume in cocktail (6rxn)	Final conc. in 1 rxn (1ul)
Sterile DW		25.75	6	154.5	Add to 50 μl
PCR buffer	10x	5	6	30	1x
dNTPs	10mM	1	6	6	200 μM
MgCl ₂	25mM	3	6	18	1.5 mM
<i>E. faecalis</i> Forward primer	10 nmol/μL (10 μM)	2.5	6	15	0.5 μM
<i>E. faecalis</i> Reverse primer	10 nmol/μL (10 μM)	2.5	6	15	0.5 μM
<i>E. faecium</i> Forward primer	10 nmol/μL (10 μM)	2.5	6	15	0.5 μM
<i>E. faecium</i> Reverse primer	10 nmol/μL (10 μM)	2.5	6	15	0.5 μM
Taq polymerase	5 unit/μL	0.25	6	1.5	1.25 units
DNA template	25 ng/μL	5	-	-	50 ng
Total	-	50	-	45x6=270	-

2

Campylobacter species

Component	Stock conc.	Volume in 1 rxn (µL)	No. of rxn	Volume in cocktail (6rxn)	Final conc. In 1 rxn
Sterile DW		18.75	6	112.5	Add to 50 µl
PCR buffer	10x	5	6	30	1x
dNTPs	10mM	1	6	6	200 µM of each
MgCl ₂	25mM	3	6	18	1.5 mM
<i>C. jejuni</i> Forward primer	10 nmol/µL (10 µM)	2.5	6	15	0.5 µM
<i>C. jejuni</i> Reverse primer	10 nmol/µL (10 µM)	2.5	6	15	0.5 µM
<i>C. coli</i> Forward primer	10 nmol/µL (10 µM)	5	6	30	1 µM
<i>C. coli</i> Reverse primer	10 nmol/µL (10 µM)	5	6	30	1 µM
23SrRNA Forward primer	10 nmol/µL (10 µM)	1	6	6	0.2 µM
23SrRNA Reverse primer	10 nmol/µL (10 µM)	1	6	6	0.2 µM
Taq polymerase	5 unit/µL	0.25	6	1.5	1.25 units
DNA template	25 ng/µL	5	-	-	125 ng
Total	-	50 µL	-	45 x6 = 270µL	-



Flow of a conventional PCR experiment




1

Make a plan



2

Prepare PCR components




3

Set PCR reaction



4

Run thermocycler



5

Run electrophoresis



6

PCR product visualization

1 PCR condition for *Enterococcus*

PCR step	Temp	Time
1. Pre-denaturation	95°C	4 minutes
2. Denaturation	95°C	30 secs
3. Annealing	55°C	1 minute.
4. Extension	72°	1 minute
Repeat step 2 to 4 for 30 cycles		
5. Final extension	72°	7 minute

2 PCR condition for *Campylobacter*

PCR step	Temp	Time
1. Pre-denaturation	95°C	6 minutes
2. Denaturation	95°C	30 sec
3. Annealing	59°C	30 sec
4. Extension	72°	30 sec
Repeat step 2 to 4 for 24-30 cycles		
5. Final extension	72°	7 minute

Set PCR cycle and run PCR



Flow of a conventional PCR experiment



1

Make a plan



2

Prepare PCR components



3

Set PCR reaction



4

Run thermocycler



5


Run electrophoresis



6

PCR product visualization

Preparation of agarose gel



- Weigh out the appropriate amount of agarose needed

- Pour the agarose into a flask
- Fill a graduated cylinder with the appropriate volume of TAE buffer.
- Pour the solution into a small flask

- Heat the solution to dissolve agarose for intervals of 15-20 seconds in a microwave oven.
- Gently swirl to disperse the contents after each interval



- Pour agarose gel solutions into the casting tray.
- Wait until the gel is set.
- Do not move the casting platform until the gel sets.

Preparing gel



Running PCR products on agarose gel

- Mix DNA with loading dye
- Load to a gel.
- Apply electricity(Run DNA on gel)
- Staining
- Visualize DNA
- PCR products are visible as different “bands”.





Flow of a conventional PCR experiment



1

Make a plan



2

Prepare PCR components




3

Set PCR reaction



4

Run thermocycler



5


Run electrophoresis




6

PCR product visualization

Visualizing PCR products



UV transilluminator




DNA gel

Visualizing PCR products




Visualizing PCR products



PCR results from identification of *Campylobacter* spp.

Visualizing PCR products



PCR results from identification of *Enterococcus* spp.

Thank you



International
Vaccine
Institute



กรมวิทยาศาสตร์การแพทย์
Department of Medical Sciences

This programme is being funded by the UK Department of Health and Social Care.
The views expressed do not necessarily reflect the UK Government's official policies.